

*Identify, diagnose
and treat balance disorders*

Why use our VSR balance assessment platform?

Create a clear pathway to medical necessity for VNG studies utilizing the VSR™ by Neurocom. Achieve verifiable, consistent results while removing all subjectivity from your patients balance assessments.

Evaluate each phase of the vestibular system with the Neurocom VSR™ platform

- ✓ Eyes open on firm surface – Tests the visual, vestibular and proprioceptive inputs
- ✓ Eyes closed on a firm surface – Tests the vestibular and proprioceptive inputs
- ✓ Eyes open on a perturbed (foam) surface – Tests the visual and vestibular inputs
- ✓ Eyes closed on a perturbed (foam) surface – Tests the vestibular inputs only

Evaluate every patient over 55 as part of a comprehensive falls prevention strategy.

- ✓ Quick and easy test that takes just 5-10 minutes to complete
- ✓ Excellent tool to evaluate patients fall risk
- ✓ Results are given as normal or abnormal
 - Normal – No further testing is warranted for the current year
 - Abnormal – The patient completes our 20-question survey to eliminate any cause(s) for an abnormal result like medications, vision or hearing problems, ankle, knee or hip problems plus many other potential complicating factors
- ✓ If no complicating or contributing factors are identified by the physician then a VNG test is warranted and the prescription is written.

Secure Health, Inc.,

1729 Red Oak Run, Fort Wayne, IN 46804

Office: 260-436-1436, Cell: 260-804-4041 Fax: 928-543-4324, DanielScherer@comcast.net

vngtechnology.com